

# K2 Feather-Lite™ Foot

## Success In Every Step



The K2 Feather-Lite™ Foot combines total safety with natural, ultra-lite ease of gait for K2 patients

- ✓ 3 Resistance Settings for Dorsi and Plantar Flexion
- ✓ 2 Keel Stiffness Options: Firm or Natural
- ✓ Feather-Lite (205g to 429g)

SPECS	Pyramid Ankle	Multi-Axis Ankle	Multi-Axis HW Ankle
Weight Limit	300lb (136kg)	220lb (100kg)	300lb (136kg)
Weight	7.3oz (205g)	13.9oz (394g)	15.2oz (429g)
Sizes	22cm - 28cm R / L	22cm - 28cm R / L	22cm - 28cm R / L
Build Height	1 <sup>15</sup> / <sub>16</sub> in (5cm)	3 <sup>3</sup> / <sub>4</sub> in (9.4cm)	3 <sup>3</sup> / <sub>4</sub> in (9.4cm)
Heel Height	3/ <sub>8</sub> in (10mm)	3/ <sub>8</sub> in (10mm)	3/ <sub>8</sub> in (10mm)

Suggested L-Codes\* (Multi-Axis Ankle):

**L5972** **L5986**

Suggested L-Code\* (Pyramid Ankle):

**L5972**

### Adjustable Multi-Axis Ankle or Pyramid Ankle



Multi-Axis Ankle  
w/ Natural Keel

Pyramid Ankle  
w/ Firm Keel

### Firm or Natural Roll-Over

### 2 Foot Shell Color Options



Caucasian



Light Black

### ANKLE ADJUSTMENTS




DAW®  
Industries

Orders (800) 252-2828 • Technical Support (800) 242-8669 • [www.daw-usa.com](http://www.daw-usa.com)

© Copyright 2021, DAW Industries, San Diego, CA - All Rights Reserved


# Adjusting Plantar & Dorsi-Flexion

## of Your K2 Feather-Lite™ Foot w/ Multi-Axial Ankle



**NOTE:** Foot comes preassembled. Medium Loctite was applied to the Foot Nut. The Foot Nut was then properly torqued to 9ft/lbs (12.2 Nm). Should you wish to replace the bumpers, it is suggested that you heat up the Foot Nut to allow for its removal. Re-assembly requires medium Loctite and proper torque.


1.



Remove the Foot Nut from the multi-axis ankle

**DO NOT USE VISE GRIPS.**

2.



Replace Bumpers

**3 Different Bumper Durometers:**

**Green** = Soft    **Yellow** = Medium    **Red** = Hard


3.



Reattach Ankle to Keel, Apply medium Loctite to bolt threads.

Be sure to **Torque Foot Nut to 9 ft/lb (12.2 Nm)**

## Pyramid Ankle



Apply medium Loctite.

**Torque Foot Bolt to 9 ft/lb (12.2 Nm)**