

DAWFLO™ SOCK



CLAIM TO FAME

- Circulation enhancement through gradient pressure
- Helps venous return, capillary stimulation when properly fitted
- Helps keep skin dry
- Frictionless ensemble without piston action when used with the **DAWFLO** Sheath
- Unique limb/prosthesis bond
- Can be machine washed and tumble dried
- No shrinking, retains softness, elasticity and gradient pressure

Worldly acclaimed as the finest sock, it provides the ultimate care of the residual limb by helping promote venous return. A technical breakthrough together with old-fashioned workmanship in hand-sewing brings a matchless level of quality. The **DAWFLO** Sock is one of the greatest advancements in prosthetic care – it is not “just another sock”.

The **DAWFLO** Sock has a soft side, made of pure wool, that “clings” to the outside of the **DAWFLO** Sheath. The **DAWFLO** Sheath and Sock become integrated, yet keep their individual functions. An outer surface, made of synthetic fibers, provides better adherence to the inside of the socket.

When properly fitted, it supports the skin and muscle tissue of the residual limb. The **DAWFLO** Sock is highly effective in helping to create a perfect, total fit of the prosthesis to the limb. Because the inside is made of the finest 100% virgin wool, the **DAWFLO** Sock is very efficient in absorbing perspiration. Available in 3, 5 and 6 ply. The special design allows for easy circulation of the air, thus evaporating perspiration. Besides wool, it can also be ordered in Cotton.

HOW TO ORDER

Please refer to the chart in Section 1.1F, p.02



DAW
INDUSTRIES