

every time

SLIP INTO YOUR SUCTION SOCKET

SMOOTHLY, NATURALLY



(No Talcum, Cream or Liquid Necessary)

- Avoid tissue misplacement, the common curse of traditional donning methods.
- Eliminate tissue damage and skin irritation.
- Ensure ideal donning and perfect positioning in the suction socket every time.

The ideal solution for the Above-Knee Suction Socket Wearer. Any other means of donning a properly fitted AK socket encourages misplacement of tissues.

The **E.Z. IN's** material is so slippery it will convince every patient that donning a tight AK Socket is not a big deal!



HOW TO ORDER

Stock #: **GEZIN-x**







x = size (Please Refer to Chart on Back)

Each E.Z. IN comes in a convenient pocket pouch.



DAW INDUSTRIES





HOW TO ORDER

	ISCHIAL CIRCUMFERENCE	Under 19 in (48 cm)	19 to 23 in (48 to 58 cm)	23 to 27 in (58 to 69 cm)	27 to 31 in (69 to 79 cm)
	SIZE	Small	Medium	Large	Extra Large
	STOCK #	GEZIN-S	GEZIN-M	GEZIN-L	GEZIN-XL



DAW INDUSTRIES