

Cool Liner™

Patented

♥ Love
the
difference!®

Patented
Vertical
Stretch
Control

Also available
without V.S.C.

**NEW
HIGH
ACTIVITY**

**ABSOLUTELY !
NO STRETCH •**
(VERTICALLY)

**THE ULTIMATE
HEAT FORMABLE GEL
for Socket Comfort**

- ✓ Maximum Cushion Comfort
- ✓ Total Sheer Protection
- ✓ Guaranteed Durability
- ✓ Mineral Oil Based Gel
- ✓ Vitamin A, D & E Enriched



**DAW®
INDUSTRIES**

© Copyright 2015, DAW Industries, San Diego, CA. All Rights Reserved.

Orders (800) 252-2828 • Technical Hotline (877) 242-2423 • www.daw-usa.com




Cool Liner™

The **ULTIMATE** Socket Comfort



DISTAL CIRCUMFERENCE



16/19
Small Uniform

18/22
Medium Uniform

20/28
Med+ Uniform

23/32
Large Uniform

28/34
Large+ Uniform

33/40
X-Large Uniform

- Ultimate Sheer Protection
- Vitamin A, D, & E Enriched
- Guaranteed Durability
- Mineral Oil Based Gel
- Heat Formable

1 When residual limb is extremely pointed (high MTP corc. with low distal circ.), you may use:

- 1** = 17/22 (Medium Tapered)
- 2** = 18/28 (Medium+ Tapered)

How to use this Size Chart

1. Measure the distal circumference of the patient's residual limb. Find the nearest number in the blue top row of the chart. Identify the column where that number is.
2. Measure the MTP circumference of your patient's residual limb. Find the nearest number in the **green** far left column of the chart. Identify the row where that number is.
3. The colored zone where the column and the row meet is the Cool Liner size (for example: column 28 and row 38 "meet" in colored zone "Med + Uniform". Therefore, the desired size is: "Med+Uniform").

Orders (800) 252-2828

Technical Hotline (877) 242-2423

www.daw-usa.com