

AK

DAWREDUCER™

with WAIST BELT

PERFECT AK EDEMA REDUCTION IS KIND TO THE STITCHES!



L-CODE : L8460 – Prosthetic Shrinker, AK

CLAIM TO FAME

- An efficient and reliable method of reducing residual limb edema for AK extremity amputees.
- Guarantees, day after day, the optimum edema reduction.
- Eliminates the risk of improper bandaging.
- Supercedes and makes obsolete all other “stump shrinker” methods.
- Easily adjustable Suspension Belt included to prevent slippage.
- Very thin, yet very strong elastic material allows for custom modification.



The **DAWREDUCER** is easy to use even for arthritic hands. The contour-weave prevents medial chafing and abrasion problems. It always provides the right compression when properly fitted. Its special lining prevents downward slippage.

HOW TO ORDER

Stock #: **SSS-1AK-x**

x = size (Please Refer to Chart on Back)



**DAW
INDUSTRIES**

HOW TO DETERMINE THE RIGHT SIZE:


AK

DAWREDUCER™ SIZING



ABOVE KNEE

Take the circumference of the thigh (in inches) 2" below the ischial level – Find corresponding size on chart below:

	SIZE	2	3	4	5	6
 <p>thigh</p>	<p>AK THIGH circumference 2" below ischial level</p>	15 in	17 in	19 in	21 in	24 in

HIGHLY RECOMMENDED
To Apply, Don't Forget to Use the DAW'NER



See DAW'NER in Section 1.4B p.01



DAW INDUSTRIES

