

ADULT

1

MULTI-AXIS ADD-ON ANKLES

Can be used with any Rhythm Foot™ or SACH Foot

For additional information on Rhythm Foot, see Section 7A, p.01.



Nut ØØ2



Nut ØØ2

Stock Number							
TAJ-MPYRNUT _{XXX} *	AA			11 1/8 oz (318 g)	2 5/8 in (6.2 cm)	See chart below	A B
TAJ-MPYRØØ2 _{XXX} *				10 7/8 oz (312 g)			

* Patients up to 130 lbs (60 kg) xxx = Ø60
 " " 175 lbs (80 kg) xxx = Ø80
 " " 220 lbs (100 kg) xxx = 1Ø0

Stock Number							
TAJ-CLNUT _{XXX} *	AA			12 3/8 oz (356 g)	3 1/16 in (81 mm)	See chart below	B
TAJ-CLØØ2 _{XXX} *				12 oz (346 g)			

2

GRAPH-LITE™ MULTI-AXIS ANKLES

Must be used with TFF-Ø1_{XXR/L}

For additional information on TFF-Ø1, see Section 7B, p.02.



Stock Number							
TGA-ØM (Adjustable Heel Height)	G			7 1/4 oz (232 g)	3 1/8 in (78 mm)		A B

SYMBOL LEGEND™

Material		Attachment Method		Dimensions		Maximum Patient Weight	
A = Aluminum Alloy	P = Plastic			H = Ht.	W = Wt.	Ø = Dia.	80 lb (36 kg)
AA = Aerospace Aluminum	S = Stainless						220 lb (100 kg)
F = Foam	T = Titanium						120 lb (54 kg)
G = Graph-Lite™	W = Wood						165 lb (75 kg)
							243 lb (110 kg)
							250 lb (113 kg)
							200 lb (90 kg) high activity



DAW INDUSTRIES

