


# FINISHINGS



## SECTION 9.2

SECTION	PAGE(S)
<b>9.2A.1 DAWSkin™</b>	
• BK DAWSKIN, BK DAWToeSKIN . . . . .	01 - 04
• AK 1- Part, DAWSKIN. . . . .	05
• AK Discontinuous 2-Part, DAWSKIN . . . . .	06
• DAWSKIN Accessories / Options	
• S-Glue . . . . .	09
• DAWSKIN Cleaner . . . . .	09
• Rubber Bands . . . . .	10
• DAWTrim™. . . . .	10
• Custom Hair . . . . .	10
<b>9.2A.2 Shower Seal™</b>	
• Shower Seal™, anti-slip, waterproof . . . . .	11
• The Story Behind this amazing innovation . . . . .	12
• Skin Tone selection. . . . .	13
• Ordering the right size. . . . .	14
<b>9.2B Cosmetic Stockings</b> . . . . .	01



# DAWSKIN™

## #1 IN THE WORLD!

**NEW**  
**MEGA-STRETCH™**

Applies like a stocking,  
lasts like a DAWSkin!



**NEW**  
**ANKLE**  
**STRETCH**  
FOR ENERGY FEET

### Off The Shelf

- Removable, flexible, and durable
- Skin-like texture and feel
- Memory-STRETCH® to easily don
- Medicare reimbursed



### 10 Skin Tones




**DAW**  
**INDUSTRIES**



# HOW TO ORDER THE RIGHT DAWSKIN

<b>DAWSKIN STYLE</b> <i>Required</i>	<b>TOE STYLE</b> <i>Required</i>	—	<b>AK/BK</b> <i>Required</i>	<b>SIZE</b> <i>Required</i>	<b>COLOR</b> <i>Required</i>	<b>FOOT SIDE</b> <i>Optional</i>
<input type="text"/>	<input type="text"/>		<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>



A	Item ID
ORIGINAL	DSK
MegaStretch	DSM
EZ Access	DSE

B	Toe Style
No Toes	Blank
Toes	T
Split Toe	S

C	AK/BK
	AK
	BK
	2PAK

E	Color
Blush	BLUSH
Fair	FAIR
Shell	SHELL
Cream	CREAM
Caucasian	CAUC
Suntan	SUN
Dark Suntan	DKSUN
Light Black	LTBLK
Black	BLACK
Graphite Black	GRHBLK

F	Only for Toes
Left	L
Right	R



D	AK Top	Calf	Ankle	Foot Size	Length AK	Length BK
<b>Adolescent</b>	14.25 - 16.25	9.125 - 10.25	5.75 - 7.00	14 - 19cm	27	20
<b>Petite</b>	17 - 19	11.75 - 13.00	7.25 - 8.50	20 - 25cm	29	20
<b>Small</b>	17 - 19	12.75 - 14.25	8.50 - 9.50	21 - 24cm	29	21
<b>Medium</b>	19 - 23	14.00 - 16.25	9.00 - 10.25	25 - 27cm	31	23
<b>Large</b>	21 - 27	16.50 - 18.75	9.75 - 11.50	28 - 30cm	34	24

**Part Number Example**

The Prosthetist would like to purchase a MegaStretch DAWSkin for a Petite BK, no toes, Caucasian Color

The Prosthetist would use the following part number format to submit the order

<b>DAWSKIN STYLE</b>	<b>TOE STYLE</b>	—	<b>AK/BK</b>	<b>SIZE</b>	<b>COLOR</b>	<b>FOOT SIDE</b>
DSM	(BLANK)		BK	P	CAUC	(BLANK)



**DAW INDUSTRIES**

# PROTECT YOUR ART

## with Lasting Beauty

Need to go back in?  
Do it with DAWSkin  
See Section 9.2A

Now with or without Toes  
or with Split Toes!

### Off the Shelf

## DAWSKIN™ or DAW<sup>toe</sup>SKIN™

or *Split toe*



#### L-CODES:

L 5962 - BK Addition, Flexible Protective Outer Surface Covering System

L 5964 - AK Addition, Flexible Protective Outer Surface Covering System

- 9 Skin Tones + Graphite Black
- 5 BK Sizes
- 5 Continuous AK Sizes
- 3 Discontinuous 2-Part AK Sizes



DAW<sup>toe</sup>SKIN

DAW<sup>toe</sup>SKIN with Hair

DAWSKIN



**DAW**  
INDUSTRIES

**DAWSKIN™**

**EZ access™**  
**DAWSKIN™**

**DAW<sup>⚡</sup>toe SKIN™**

**DAW<sup>Split</sup>toe SKIN™**

**BK**

# Off the Shelf

The Custom Flexible Protective Skin System Claim to fame



- Slips on and off easily
- Soft, flexible and durable
- Skin-like texture and feel
- 3 Toe-Styles fit most manufacturers' feet
- Optional custom Hair for DAWSkin, DAWtoeSKIN & DAW<sup>Split</sup>toe SKIN

9 Skin Tones + 1 Graphite Black

5 BK DAWSkin Sizes

Adolescent, Petite, Small, Medium and Large

4 BK DAWtoe Skin & DAW<sup>Split</sup>toe Skin Sizes

Petite, Small, Medium and Large

**L-CODE:**

L 5962 – BK Addition, Flexible Protective Outer Surface Covering System



**Sizing Definition**

Adolescent	(14-19)*
Petite	(20-25)**
Small	(21-24)*
Medium	(25-27)*
Large	(28-30)*

\* Foot sizes in centimeters (See page 9.2A p.02)  
\*\* Designed for the slender foot and ankle



**DAW INDUSTRIES** See Color Swatch on page Section 9.2A.1 p.03



# EZ access™ DAWSKIN™

“PROTECTIVE OUTER SEAL” ©

Section 9.2A.1 p.05

**NEW**

**EZ WRAP™**  
matching DawSkin material

2-part  
included  
in box

PATENTED



*On & Off  
like a sock!*

**EZ access™  
DAWSkin**

## Off The Shelf

- Removable, flexible, and durable
- Skin-like texture and feel
- Memory-STRETCH® for easy donning
- Medicare reimbursed L 5962

**NEW  
TEAR-FREE  
FLEX  
ANKLE**







**10 Skin Tones**



**DAW®  
INDUSTRIES**

# DAWSKIN & DAW<sup>toe</sup>SKIN STANDARD & MEGA-STRETCH

## SIZING CHART



- BK, AK and 2-Part AK to choose from.
- Slips on and off easily.
- Soft, flexible and durable.
- Skin-like texture and feel.
- Does not interfere with function of knee.
- **Nine skin tones + Graphite Black:**  
(Blush, Fair, Shell, Cream, Caucasian, Suntan, Dark Suntan, Light Black, Black, Graphite Black)

SKIN SIZE CIRC.	ADOLESCENT	PETITE	SMALL	MEDIUM	LARGE
	14 to 19 cm	20 to 25 cm	21 to 24 cm	25 to 27 cm	28 to 30 cm
<b>A</b>	14 1/8 to 16 1/4 in (35 to 41 cm)	17 to 19 in (43 to 48 cm)	17 to 19 in (43 to 48 cm)	19 to 23 in (48 to 58 cm)	21 to 27 in (52 to 68 cm)
<b>B</b>	9 1/8 to 10 1/4 in (23 to 26 cm)	11 3/4 to 13 in (29 to 33 cm)	12 3/4 to 14 1/4 in (32 to 36 cm)	14 to 16 1/4 in (35 to 41 cm)	16 1/2 to 18 3/4 in (41 to 47 cm)
<b>C</b>	5 1/4 to 7 in (13 to 18 cm)	7 1/4 to 8 1/2 in (18 to 21 cm)	8 1/2 to 9 1/2 in (21 to 24 cm)	9 to 10 1/4 in (22 to 26 cm)	9 3/4 to 11 1/2 in (24 to 29 cm)
<b>LENGTH</b>	<b>BK</b> 20 in (50 cm) <b>AK</b> 27 in (68 cm)	<b>BK</b> 20 in (50 cm) <b>AK</b> 29 in (73 cm)	<b>BK</b> 21 in (53 cm) <b>AK</b> 29 in (73 cm)	<b>BK</b> 23 in (58 cm) <b>AK</b> 31 in (78 cm)	<b>BK</b> 24 in (60 cm) <b>AK</b> 34 in (85 cm)

PRODUCT AVAILABILITY	BK	AK	AK2-PART
DAWSKIN™ ORIGINAL	ALL SIZES	N/A	N/A
DAWSKIN MEGASTRETCH™	ALL SIZES	ALL SIZES	ALL SIZES
DAW <sup>toe</sup> SKIN™ MEGASTRETCH	ALL SIZES (except adolescent)	ALL SIZES (except adolescent)	ALL SIZES (except adolescent)
DAW <sup>splittoe</sup> SKIN™ MEGASTRETCH	ALL SIZES (except adolescent)	ALL SIZES (except adolescent)	ALL SIZES (except adolescent)
DAWSKIN EZ-access™ MEGASTRETCH <small>-Toes -Split toes -No toes</small>	ALL SIZES	N/A	N/A



**DAW INDUSTRIES**



# DAWSKIN™

## AK 1-PART

**Off the Shelf**  
The Custom Flexible Protective Skin System



Claim to fame



- Slips on and off easily
- Soft, flexible and durable
- Skin-like texture and feel
- Does not interfere with function of knee

**9 Skin Tones + Graphite Black**

Blush • Fair • Shell • Cream • Caucasian  
Suntan • Dark Suntan • Light Black • Black

**5 AK Sizes**

Adolescent (14-19)\*, Petite ° (22-25)\*, Small (21-24)\*,  
Medium (25-27)\*, and Large (28-30)\*.

\* Foot sizes in centimeters (See page 9.2A p.02)  
° Designed for the slender foot and ankle

**L-CODE:**

L 5964 – AK Addition, Flexible Protective  
Outer Surface Covering System



**HOW TO ORDER**

Stock #:  
**DSM-AK<sub>xy</sub>**

x = size      y = skin tone

See Color Swatch on page Section 9.2A.1 p.09



**DAW INDUSTRIES**







# DAWSKIN™

## AK DISCONTINUOUS 2-PART



**Off the Shelf**  
 The Custom Flexible Protective Skin System



**Claim to fame**

- Superior Durability
- Ultra-Light
- Enhances sitting cosmesis
- Eliminates foam cover stress during flexion
- No more torn or ripped foam covers
- Pre-drilled to accommodate knee unit and 30mm pylon
- Fits all Stable-Lite 4 bar knees (refer to section B – “AK System”)

**9 Skin Tones + Graphite Black**

Blush • Fair • Shell • Cream • Caucasian  
 Suntan • Dark Suntan • Light Black • Black

**4 AK Sizes**

Petite (20-25) Small (21-24)\*, Medium (25-27)\*,  
 Large (28-30)\*

*\*Foot sizes in centimeters (See page 9.2A p.00b)*



Easy access to back of AK for knee adjustments.

**SYSTEM INCLUDES:**

- Shin Section
- Femoral Section
- Valve Hole Reinforcement

**L-CODE:**

L 5964 – AK Addition, Flexible Protective Outer Surface  
 Covering System

**HOW TO ORDER**

Stock #:  
**MegaStretch: DSM-2PAKxy**  
 Also **DSMT-2PAKxy**  
**DSMS-2PAKxy**



**DAW INDUSTRIES**





# DAWSKIN™ ACCESSORIES



## S-GLUE™ : REGULAR FORMULA



The only sure, long lasting tough glue that holds the **DAWSKIN** to your socket without altering the quality of the **DAWSKIN**. Can also be used as a fast setting glue for 1000 other applications.

### HOW TO ORDER

Stock #:  
DSK-SG

## S-GLUE™ : RELEASE FORMULA



For gluing **DAWSKIN** to any material (sockets, plastics, rubber, wood, leather, foam, etc.). Ideal for repairs.

### HOW TO ORDER

Stock #:  
DSK-SGREL

## Socket Trim




25 meter roll of our most popular socket edging material. Uncoated for long term durability.

### HOW TO ORDER

Stock #:  
FP-ST: 25 Meter Roll

## DAWTRIM™



The most elegant finish for the socket brim. **DAWTRIM** is a highly elastic strip specially designed for custom finishing the socket brim with a cushioned edge. Available now in all ten **DAWSKIN** colors.

### HOW TO ORDER

Stock #:  
DSK-DT: Enough for 1 socket  
DSK-DT5YD: 5 yards

**9 Skin Colors Available + Graphite Black**  
Blush, Fair, Caucasian, Cream, Shell, Suntan, Dark Suntan, Light Black, Black, Graphite Black

*(specify color)*

## DAWSKIN CLEANER



**DAWSKIN** can be cleaned with soap and water – Some hard to remove stains need the **DAWSKIN** cleaner. Sold individually.

### HOW TO ORDER

Stock #:  
DSK-CLEANER

## Color Swatch



Years of experience prove this selection of 9 skin tones + Graphite Black agrees with all patients skin tones. More choices prove more confusing than helpful.

### HOW TO ORDER

Stock #:  
DSK-CS

*The skin tones are:* Blush, Fair, Shell, Cream, Caucasian, Suntan, Dark Suntan, Light Black, Black, Graphite Black




**DAWSKIN™**

**OPTIONS**




**RUBBER BANDS**




When frequent removal of the **DAWSKIN** is expected, one might choose to glue **DAWSKIN** to a pre-installed DAW proximal Rubber Band rather than directly to the socket. To do so, select correct Rubber Band by measuring circumference of prosthesis below Patella trim line:

**HOW TO ORDER**

Circumference	12 - 14 in	14 - 17 in	17 - 20 in	20 - 23 in
Stock #  (Package of 1)	TRB-ØB	TRB-ØAM	TRB-ØAL	TRB-ØAX

**SUCTION VALVE RUBBER CUP**




The clean way to finish the suction valve hole in the foam cover. Use S-Glue (DSK-SG) to install.

**HOW TO ORDER**

Stock #:  
TRC-Ø

**VALVE PATCH**



A perfect complement to the Suction Valve Rubber Cup. For a total finish, glue DAW's Valve Patch with S-Glue to the AK portion of the foam covering.

**HOW TO ORDER**

Stock #:  
DSK-VP

**Custom HAIR**



Most men prefer the natural look of hair. Any **DAWSKIN** can be custom finished this way. This is not "off the shelf" as all our other **DAWSKIN** products, therefore allow 72 hours for shipping.

**HOW TO ORDER**

DSM-CUSTOMHAIR

(One hair color only, other colors do not show.)



**Be Safe!**  
Just like before  
**Shower**  
with  
confidence on  
your own 2 legs!

Patent  
Pending




# Shower Seal™

“Anti-slip”  
Waterproof

Be it...

- in the shower
- bathing
- by the pool

Enjoy the safety and  
ease of standing on  
both your legs.



Link to DAW Shower Seal Video  
<https://goo.gl/1mxcUH>



**DAW**  
**INDUSTRIES**

# THE STORY BEHIND THIS AMAZING INNOVATION

Dear Practitioner,

I remember in my childhood seeing my father David, a World War II BK amputee, hopping around on one leg to take a shower in our family bathroom. I also remember hearing on occasion a thundering noise... the scary sign of my father slipping and falling down! Thank God he was never seriously hurt.


Then I thought it was just a normal occurrence. No more! With DAW's waterproof anti-slip Shower Seal, every amputee can now navigate on both feet as safely as a non-amputee in any wet environment (bathroom, swimming pool, beach, etc.).

Our Patent Pending Daw Shower Seal comes in many skin tones to be less conspicuous in any surroundings. Here is a convenient selection swatch for the Shower Seal which will help all your BK amputees live a safer daily life and feel whole again.



**DAW INDUSTRIES**






# Shower Seal™

## “Anti-slip”

### Waterproof

- **Be Safe!**
- **Shower with confidence on your own 2 legs!**

The Tibial Section of the Shower Seal Matches your patient’s skin tone.



**DAW**  
**INDUSTRIES**



# HOW TO ORDER THE RIGHT TIBIAL SHOWER SEAL

SHOWER SEAL SYSTEM	STYLE <i>Required</i>	SIZE <i>Required</i>	COLOR <i>Required</i>	FOOT SIDE <i>If Toes Ordered</i>
SSS	<b>A</b>	<b>B</b>	<b>C</b>	<b>D</b>

A Style	
No Toes	Ø
Toes	T
Split Toe	S

D Only for Toes	
Left	L
Right	R

B	Foot Size		Calf	Ankle	Length *	Size Code
Adoles.	14 - 19 cm	cm	23 - 26	14.5 - 17.75	50.8	A
		in	<b>9 - 10</b>	<b>5.5 - 7</b>	<b>20</b>	
Petite	20 - 25 cm	cm	30 - 33	18.5 - 21.5	50.8	P
		in	<b>11.75 - 13</b>	<b>7.25 - 8.5</b>	<b>20</b>	
Small	21 - 24 cm	cm	32.25 - 36.25	21.5 - 24	53.34	S
		in	<b>12.75 - 14.25</b>	<b>8.5 - 9.5</b>	<b>21</b>	
Med.	25 - 27 cm	cm	35.5 - 41.25	22.75 - 26	58.42	M
		in	<b>14 - 16.25</b>	<b>9 - 10.25</b>	<b>23</b>	
Large	28 - 30 cm	cm	41.75 - 47.75	24.75 - 29.25	60.96	L
		in	<b>16.5 - 18.75</b>	<b>9.75 - 11.5</b>	<b>24</b>	

C Color	
Blush	BLUSH
Fair	FAIR
Shell	SHELL
Cream	CREAM
Caucasian	CAUC
Suntan	SUN
Dark Suntan	DKSUN
Light Black	LTBLK
Black	BLACK
Graphite Black	GRHBLK

\* Maximum length floor to top of socket

**Part Number Example**

The Prosthetist would like to purchase a Shower Seal System for a Petite BK, no toes, Caucasian Color

The Prosthetist would use the following part number format to submit the order

SHOWER SEAL SYSTEM	STYLE <i>Required</i>	SIZE <i>Required</i>	COLOR <i>Required</i>	FOOT SIDE <i>If Toes Ordered</i>
SSS	Ø	P	CAUC	



**DAW INDUSTRIES**

**Must Use: Daw S-Glue**  
 Good enough to finish up to 18 BK prosthesis!  
 (See page 4 paragraph 10)



**S-Glue DSK-SG**

# COSMETIC STOCKINGS

## COSMETIC STOCKINGS TO PROTECT AND BEAUTIFY



DAW Sheer Stocking

DAW Ribbed Stocking

### Claim to fame

- Easy to apply
- Very inexpensive
- 3 lengths AK and BK fit all
- Tubular design fits any size foot



The DAW Ribbed Stocking substantially prolongs the life of any Foam Cover. Apply the ribbed stocking(s) first. Then apply sheer stocking, which helps to achieve a more natural looking prosthesis. Since these are tube stockings, please take into account the foot length when measuring your patient.

### HOW TO ORDER

<b>TCS-SA<sub>xx</sub></b> <b>TCS-SB<sub>xx</sub></b>	<b>TCS-RA<sub>xx</sub></b> <b>TCS-RB<sub>xx</sub></b>
--	--

**Lengths available**

**A** = AK..... xx = 35, 38, and 41 in (89, 95, and 104 cm)

**B** = BK..... xx = 18, 21 and 24 in (45, 53, and 60 cm)

xx = Length



**DAW**  
**INDUSTRIES**





## Remember to look for...



This seal of excellence is earned only by very select items. Although all **DAW** products meet the highest standard of quality, the “Signature of Excellence” items address themselves to the select patients who demand the best of the best.



**DAW**  
**INDUSTRIES**

© Copyright 2015, DAW Industries, San Diego, CA. All Rights Reserved.

Orders (800) 252-2828 • Technical Hotline (877) 242-2423 • [www.daw-usa.com](http://www.daw-usa.com)

