

Be Safe!
Just like before
Shower
with
confidence on
your own 2 legs!


**Patent
Pending**



Shower Seal™

"Anti-slip"

Waterproof



Link to DAW Shower Seal Video
<https://goo.gl/1mxcUH>

Be it...

- in the shower
- bathing
- by the pool

**Enjoy the safety and
ease of standing on
both your legs.**



**DAW
INDUSTRIES**

© Copyright 2015, DAW Industries, San Diego, CA. All Rights Reserved.

Orders (800) 252-2828 • Technical Hotline (877) 242-2423 • www.daw-usa.com

THE STORY BEHIND THIS AMAZING INNOVATION

Dear Practitioner,

I remember in my childhood seeing my father David, a World War II BK amputee, hopping around on one leg to take a shower in our family bathroom. I also remember hearing on occasion a thundering noise... the scary sign of my father slipping and falling down! Thank God he was never seriously hurt.

Then I thought it was just a normal occurrence. No more! With DAW's waterproof anti-slip Shower Seal, every amputee can now navigate on both feet as safely as a non-amputee in any wet environment (bathroom, swimming pool, beach, etc.).

Our Patent Pending Daw Shower Seal comes in many skin tones to be less conspicuous in any surroundings. Here is a convenient selection swatch for the Shower Seal which will help all your BK amputees live a safer daily life and feel whole again.



Shower Seal Sleeve
(1)



Foam Cover
(not included)



Tibial Shower Seal
(2)



Finished Prosthesis with Shower Seal System

*** Shower Seal Box Includes: 1, 2, 3 and 4**



Latex Strap
(3)




2x Quick Fix Patches
(4)



DAW INDUSTRIES

© Copyright 2015, DAW Industries, San Diego, CA. All Rights Reserved.

Orders (800) 252-2828 • Technical Hotline (877) 242-2423 • www.daw-usa.com




Shower Seal™

"Anti-slip"

Waterproof

- **Be Safe!**
- **Shower with confidence on your own 2 legs!**

The Tibial Section of the Shower Seal Matches your patient's skin tone.



**DAW
INDUSTRIES**

© Copyright 2015, DAW Industries, San Diego, CA. All Rights Reserved.

Orders (800) 252-2828 • Technical Hotline (877) 242-2423 • www.daw-usa.com

HOW TO ORDER THE RIGHT TIBIAL SHOWER SEAL

SHOWER SEAL SYSTEM	STYLE <i>Required</i>	SIZE <i>Required</i>	COLOR <i>Required</i>	FOOT SIDE <i>If Toes Ordered</i>
SSS —	A	B	C	D

A Style	
No Toes	Ø
Toes	T
Split Toe	S

D Only for Toes	
Left	L
Right	R

B	Foot Size		Calf	Ankle	Length *	Size Code
Ado-les.	14 - 19 cm	cm	23 - 26	14.5 - 17.75	50.8	A
		in	9 - 10	5.5 - 7	20	
Petite	20 - 25 cm	cm	30 - 33	18.5 - 21.5	50.8	P
		in	11.75 - 13	7.25 - 8.5	20	
Small	21 - 24 cm	cm	32.25 - 36.25	21.5 - 24	53.34	S
		in	12.75 - 14.25	8.5 - 9.5	21	
Med.	25 - 27 cm	cm	35.5 - 41.25	22.75 - 26	58.42	M
		in	14 - 16.25	9 - 10.25	23	
Large	28 - 30 cm	cm	41.75 - 47.75	24.75 - 29.25	60.96	L
		in	16.5 - 18.75	9.75 - 11.5	24	

* Maximum length floor to top of socket

C Color	
Blush	BLUSH
Fair	FAIR
Shell	SHELL
Cream	CREAM
Caucasian	CAUC
Suntan	SUN
Dark Suntan	DKSUN
Light Black	LTBLK
Black	BLACK
Graphite Black	GRHBLK

Part Number Example

The Prosthetist would like to purchase a Shower Seal System for a Petite BK, no toes, Caucasian Color

The Prosthetist would use the following part number format to submit the order

SHOWER SEAL SYSTEM	STYLE <i>Required</i>	SIZE <i>Required</i>	COLOR <i>Required</i>	FOOT SIDE <i>If Toes Ordered</i>
SSS —	Ø	P	CAUC	

Must Use: Daw S-Glue

Good enough to finish up to 18 BK prosthesis!
(See page 4 paragraph 10)



**S-Glue
DSK-SG**



**DAW
INDUSTRIES**

© Copyright 2015, Daw Industries, San Diego, CA. All Rights Reserved.

Orders (800) 252-2828 • Technical Hotline (877) 242-2423 • www.daw-usa.com