

 TM





Applies like a stocking, lasts like a DAWSkin!







- Removable, flexible, and durable
- Skin-like texture and feel
- Memory-STRETCH[®] to easily don
- Medicare reimbursed








© Copyright 2015, DAW Industries, San Diego, CA. All Rights Reserved.

HOW TO ORDER THE RIGHT DAWSKIN



D	АК Тор	Calf	Ankle	Foot Size	Length AK	Length BK
Adolescent	14.25 - 16.25	9.125 - 10.25	5.75 - 7.00	14 - 19cm	27	20
Petite	17 - 19	11.75 - 13.00	7.25 - 8.50	20 - 25cm	29	20
Small	17 - 19	12.75 - 14.25	8.50 - 9.50	21 - 24cm	29	21
Medium	19 - 23	14.00 - 16.25	9.00 - 10.25	25 - 27cm	31	23
Large	21 - 27	16.50 - 18.75	9.75 - 11.50	28 - 30cm	34	24

Part Number Example

The Prosthetist would like to purchase a MegaStretch DAWSkin for a Petite BK, no toes, Caucasian Color

The Prosthetist would use the following part number format to submit the order

DAWSKIN STYLE	TOE STYLE		AK/BK SIZE	COLOR	FOOT SIDE	
DSM	(BLANK)	_	BK P	CAUC	(BLANK)	



PROTECT YOUR ART

with Lasting Beauty

Need to go
back in?
Do it with DAWSkin
See Section 9.2A

Now with or without Toes or with Split Toes!

Off the Shelf DAWSKIN or DAW toe SKIN



L-CODES:


L 5962 - BK Addition, Flexible Protective Outer Surface Covering System

L 5964 - AK Addition, Flexible
Protective Outer Surface
Covering System

- 9 Skin Tones + Graphite Black
- 5 BK Sizes
- 5 Continuous AK Sizes
- 3 Discontinuous 2-Part AK Sizes

















Off the Shelf

The Custom Flexible Protective Skin System



Claim to fame

- Slips on and off easily
- Soft, flexible and durable
- Skin-like texture and feel
- 3 Toe-Styles fit most manufacturers' feet
- Optional custom Hair for DAWSkin, DAWtoeSKIN & DAWtoeSKIN
- 9 Skin Tones + 1 Graphite Black
- 5 BK DAWSkin Sizes
 Adolescent, Petite, Small, Medium and Large
- 4 BK DAW toe Skin & DAW toe Skin Sizes
 Petite, Small, Medium and Large

L-CODE:

L 5962 – BK Addition, Flexible Protective Outer Surface Covering System



Sizing Definition

Adolescent (14-19)*
Petite (20-25)**
Small (21-24)*

Medium (25-27)* Large (28-30)*

izes in centimeters



See Color Swatch on page Section 9.2A.1 p.03



^{*} Foot sizes in centimeters (See page 9.2A p.02)


^{**} Designed for the slender foot and ankle


PROTECTIVE OUTER SEAL"



DAWSKIN & DAWtoe SKIN STANDARD & MEGA-STRETCH

SIZING CHART





- BK, AK and 2-Part AK to choose from.
- Slips on and off easily.
- Soft, flexible and durable.
- Skin-like texture and feel.
- Does not interfere with function of knee.
- Nine skin tones + Graphite Black: (Blush, Fair, Shell, Cream, Caucasian, Suntan, Dark Suntan, Light Black, Black, Graphite Black)

SKIN SIZE CIRC.	ADOLESCENT 14 to 19 cm	PETITE 20 to 25 cm	SMALL 21 to 24 cm	MEDIUM 25 to 27 cm	LARGE 28 to 30 cm
A	14 1/8 to 16 1/4 in (35 to 41 cm)	17 to 19 in (43 to 48 cm)	17 to 19 in (43 to 48 cm)	19 to 23 in (48 to 58 cm)	21 to 27 in (52 to 68 cm)
В	9 ½ to 10 ½ in (23 to 26 cm)	11 ³ ⁄4 to 13 in (29 to 33 cm)	12 ³ / ₄ to 14 ¹ / ₄ in (32 to 36 cm)	14 to 16 ¼ in (35 to 41 cm)	16 ½ to 18 ¾ in (41 to 47 cm)
С	5 ½ to 7 in (13 to 18 cm)	7 ½ to 8 ½ in (18 to 21 cm)	8½ to 9½ in (21 to 24 cm)	9 to 10 1/4 in (22 to 26 cm)	9 3/4 to 11 1/2 in (24 to 29 cm)
LENGTH	BK 20 in (50 cm) AK 27 in (68 cm)	BK 20 in (50 cm) AK 29 in (73 cm)	BK 21 in (53 cm) AK 29 in (73 cm)	BK 23 in (58 cm) AK 31 in (78 cm)	BK 24 in (60 cm) AK 34 in (85 cm)

PRODUCT AVAILABILITY	ВК	AK	AK2-PART
DAWSKIN [™] ORIGINAL	ALL SIZES	N/A	N/A
DAWSKIN MEGASTRETCH™	ALL SIZES	ALL SIZES	ALL SIZES
DAWtoeSKIN [™] MEGASTRETCH	ALL SIZES (except adolescent)	ALL SIZES (except adolescent)	ALL SIZES (except adolescent)
DAWsplittoeSKIN [™] MEGASTRETCH	ALL SIZES (except adolescent)	ALL SIZES (except adolescent)	ALL SIZES (except adolescent)
DAWSKIN EZ-access [™] MEGASTRETCH -Toes -Split toes -No toes	ALL SIZES	N/A	N/A







Off the Shelf The Custom Flexible Protective Skin System







Slips on and off easily

- Soft, flexible and durable
- Skin-like texture and feel
- Does not interfere with function of knee
- 9 Skin Tones + Graphite Black

Blush • Fair • Shell • Cream • Caucasian Suntan • Dark Suntan • Light Black • Black

5 AK Sizes

Adolescent (14-19)*, Petite °(22-25)*, Small (21-24)*, Medium (25-27)*, and Large (28-30)*.

- * Foot sizes in centimeters (See page 9.2A p.02)
- Designed for the slender foot and ankle

L-CODE:

L 5964 - AK Addition, Flexible Protective Outer Surface Covering System

Also Available with Split Toes:

HOW TO ORDER

Stock #:

DSM-AKxy

x = sizey = skin toneSee Color Swatch on page Section 9.2A.1 p.09





DAW*SKIN*



DISCONTINUOUS 2-PART



Off the Shelf The Custom Flexible Protective Skin System







Easy access to back of AK for knee adjustments.

Claim to fame

- Superior Durability
- Ultra-Light
- Enhances sitting cosmesis
- Eliminates foam cover stress during flexion
- No more torn or ripped foam covers
- Pre-drilled to accommodate knee unit and 30mm pylon
- Fits all Stable-Lite 4 bar knees (refer to section B - "AK System")

9 Skin Tones + Graphite Black

Blush • Fair • Shell • Cream • Caucasian Suntan • Dark Suntan • Light Black • Black

4 AK Sizes

Petite (20-25) Small (21-24)*, Medium (25-27)*, Large (28-30)*

*Foot sizes in centimeters (See page 9.2A p.00b)

SYSTEM INCLUDES:

- Shin Section
- Femoral Section
- Valve Hole Reinforcement

L-CODE:

L 5964 - AK Addition, Flexible Protective Outer Surface Covering System

HOW TO ORDER

Stock #:

MegaStretch: DSM-2PAKxy Also DSMT-2PAKxy **DSMS-2PAK**XV





DAWSK//V ACCESSORIES



S-GLUE™: REGULAR FORMULA



The only sure, long lasting tough glue that holds the **DAWSKIN** to your socket without altering the quality of the **DAWSKIN**. Can also be used as a fast setting glue for 1000 other applications.

HOW TO ORDER

Stock #: DSK-SG

S-GLUE™: RELEASE FORMULA



For gluing **DAWSKIN** to any material (sockets, plastics, rubber, wood, leather, foam, etc.). Ideal for repairs.

HOW TO ORDER

Stock #: DSK-SGREL

Socket Trim



25 meter roll of our most popular socket edging material. Uncoated for long term durability.

Only available in Caucasian.

HOW TO ORDER

Stock #: FP-ST: 25 Meter Roll

DAW*TRIM* ™



The most elegant finish for the socket brim. **DAWTRIM** is a highly elastic strip specially designed for custom finishing the socket brim with a cushioned edge. Available now in all ten **DAWSKIN** colors.

9 Skin Colors Available + Graphite Black
Blush, Fair, Caucasian, Cream, Shell, Suntan, Dark Suntan,
Light Black, Black, Graphite Black

HOW TO ORDER

Stock #:
DSK-DT: Enough for 1
socket
DSK-DT5YD: 5 yards

(specify color)

DAWSKIN CLEANER



DAWSKIN can be cleaned with soap and water – Some hard to remove stains need the **DAWSKIN** cleaner. Sold individually.

HOW TO ORDER

Stock #:
DSK-CLEANER

Color Swatch





Years of experience prove this selection of 9 skin tones + Graphite Black agrees with all patients skin tones. More choices prove more confusing than helpful.

The skin tones are: Blush, Fair, Shell, Cream, Caucasian, Suntan, Dark Suntan, Light Black, Black, Graphite Black

HOW TO ORDER

Stock #: DSK-CS


DAWSKIN



DAWSK//V OPTIONS



RUBBER BANDS



When frequent removal of the **DAWSKIN** is expected, one might choose to glue **DAWSKIN** to a pre-installed DAW proximal Rubber Band rather than directly to the socket. To do so, select correct Rubber Band by measuring circumference of prosthesis below Patella trim line:

HOW TO ORDER

Circumference	12 - 14 in	14 - 17 in	17 - 20 in	20 - 23 in
Stock # Package of 1)	TRB-ØB	TRB-ØAM	TRB-ØAL	TRB-ØAX

SUCTION VALVE RUBBER CUP




The clean way to finish the suction valve hole in the foam cover. Use S-Glue (DSK-SG) to install.

HOW TO ORDER

Stock #: TRC-Ø

VALVE PATCH



A perfect complement to the Suction Valve Rubber Cup. For a total finish, glue DAW's Valve Patch with S-Glue to the AK portion of the foam covering.

HOW TO ORDER

Stock #: DSK-VP

Custom HAIR



Most men prefer the natural look of hair. Any **DAWSKIN** can be custom finished this way. This is not "off the shelf" as all our other **DAWSKIN** products, therefore allow 72 hours for shipping.

(One hair color only, other colors do not show.)

HOW TO ORDER

DSM-CUSTOMHAIR

