

DAWFLQ™

Section 1.1A p.01

The Choice of Professionals





DAWFLO™ ENSEMBLE

CLAIM TO FAME

- Circulation enhancement through gradient pressure
- Helps venous return, capillary stimulation when properly fitted
- Reduces or eliminates irritation and bruising
- Moisture barrier, substantially decreases perspiration
- Helps keep skin dry
- Frictionless ensemble without piston action
- Unique limb/prostheses bond
- Can be machine washed and tumble dried: No shrinking, retains softness, elasticity and gradient pressure

DAWFLO ENSEMBLE

= DAWFLO Sheath + Sock

Worldly acclaimed as the finest Sock / Sheath system, it provides the ultimate care of the residual limb by helping promote venous return. A technical breakthrough together with old-fashioned workmanship in handsewing brings a matchless level of quality. The **DAWFLO** Sock / Sheath ensemble is one of the greatest advancements in prosthetic care – it is not “just another sock”.

The **DAWFLO** Sheath has a soft surface next to the skin made of the finest sheer material. An outside surface that grips the **DAWFLO** Sock, preventing movement between Sheath and Sock. It also provides a ventilated mesh which allows perspiration to pass through.

The **DAWFLO** Sock has a soft side, made of pure wool, that “clings” to the outside of the **DAWFLO** Sheath. The **DAWFLO** Sheath and Sock become integrated, yet keep their individual functions. An outer surface, made of synthetic fibers, provides better adherence to the inside of the socket.

When properly fitted, it supports the skin and muscle tissue of the residual limb. The **DAWFLO** Sock is highly effective in helping to create a perfect, total fit of the prosthesis to the limb. Because the inside is made of the finest 100% virgin wool, the **DAWFLO** Sock is very efficient in absorbing perspiration. Available in 3, 5 and 6-ply. The special design allows for easy circulation of the air, thus evaporating perspiration. Besides wool, it can also be ordered in Cotton.

HOW TO ORDER

Please refer to the chart in Section 1.1F, p.02



**DAW®
INDUSTRIES**



DAWFLO® SHEATH™



CLAIM TO FAME

- Helps venous return, capillary stimulation when properly fitted.
- Reduces or eliminates irritation, friction and bruising.
- Moisture barrier, substantially decreases perspiration.
- Helps keep skin dry.
- Circulation enhancement through gradient pressure.
- Outlives regular Prosthetic Sheaths 10 to 1.

Worldly acclaimed as the finest Prosthetic Sheath system, it provides the ultimate care of the residual limb by helping promote venous return. A technical breakthrough together with old-fashioned workmanship in handsewing brings a matchless level of quality. The **DAWFLO** Sheath is one of the greatest advancements in prosthetic care - it is not "just another prosthetic sheath".

The **DAWFLO** Sheath has a soft surface next to the skin made of the finest sheer material. Its outside surface grips the **DAWFLO** Sock, preventing movement between Sheath and Sock. The perspiration passes through to the **DAWFLO** Sock thanks to the **DAWFLO** Sheath ventilated mesh.

When properly fitted, it supports the skin and muscle tissue of the residual limb. The **DAWFLO** Sheath is highly effective in helping to create a perfect, total fit of the prosthesis to the limb. The special design allows for easy circulation of the air, thus evaporating the perspiration and creating a moisture barrier screen which helps keep the skin dry.

HOW TO ORDER

Please refer to the chart in Section 1.1F, p.02



**DAW®
INDUSTRIES**

DAWFLO® SOCK



CLAIM TO FAME

- Circulation enhancement through gradient pressure
- Helps venous return, capillary stimulation when properly fitted
- Helps keep skin dry
- Frictionless ensemble without piston action when used with the **DAWFLO** Sheath
- Unique limb/prostheses bond
- Can be machine washed and tumble dried
- No shrinking, retains softness, elasticity and gradient pressure

Worldly acclaimed as the finest sock, it provides the ultimate care of the residual limb by helping promote venous return. A technical breakthrough together with old-fashioned workmanship in hand-sewing brings a matchless level of quality. The **DAWFLO** Sock is one of the greatest advancements in prosthetic care – it is not “just another sock”.

The **DAWFLO** Sock has a soft side, made of pure wool, that “clings” to the outside of the **DAWFLO** Sheath. The **DAWFLO** Sheath and Sock become integrated, yet keep their individual functions. An outer surface, made of synthetic fibers, provides better adherence to the inside of the socket.

When properly fitted, it supports the skin and muscle tissue of the residual limb. The **DAWFLO** Sock is highly effective in helping to create a perfect, total fit of the prosthesis to the limb. Because the inside is made of the finest 100% virgin wool, the **DAWFLO** Sock is very efficient in absorbing perspiration. Available in 3, 5 and 6 ply. The special design allows for easy circulation of the air, thus evaporating perspiration. Besides wool, it can also be ordered in Cotton.

HOW TO ORDER

Please refer to the chart in Section 1.1F, p.02



**DAW®
INDUSTRIES**

© Copyright 2015, DAW Industries, San Diego, CA. All Rights Reserved.

Orders (800) 252-2828 • Technical Hotline (877) 242-2423 • www.daw-usa.com



DAWFLO® ENSEMBLE™

KNEE DISARTICULATION



CLAIM TO FAME

- Circulation enhancement through gradient pressure.
- Helps venous return.
- Capillary stimulation.
- Moisture barrier.
- Helps skin dry.
- Frictionless ensemble without piston action.
- Unique limb/prostheses bond.
- Available in size 2-3-4-5 from 18" to 22" long.

Completing the family of DAWFLO® Socks & DAWFLO® Sheaths

The knee disarticulation amputee has always had special needs and requirements due to the shape and length of their residual limb. This contoured Sock and Sheath provides all the benefits of the DAWFLO® combination system with the added advantage of a perfect fit without the need of custom factory ordering.

HOW TO ORDER

Please refer to the chart in Section 1.1F, p.02



**DAW®
INDUSTRIES**



ENSEMBLE WITH V-BAND GARTER



ENSEMBLE WITH GARTER ONLY



DAWFLOTM

with V-BAND GARTER

or with GARTER only

CLAIM TO FAME

- Total prevention of AK sock slippage in the socket.
- Medial & Lateral attachment.
- Helps prevent skin irritation medially due to sock slippage.
- Enhances AK Residual Limb suspensory.
- All the benefits of the **DAWFLO** Ensemble.

Designed especially for the Hip-Joint prosthesis, the **DAWFLO** with V-Band Garter provides unexcelled comfort. The exclusive design reduces bulk, aids suspension, and eliminates slippage in the socket, both medially and laterally.

Patients not needing the V-band can also reduce sheath/sock slippage by using the **DAWFLO** Garter. In both cases the ensemble is kept in the pulled up position by an elastic band, folding the top over the socket brim and tying the draw string.

HOW TO ORDER

Please refer to the chart in Section 1.1F, p.11



DAW[®]
INDUSTRIES


© Copyright 2015, DAW Industries, San Diego, CA. All Rights Reserved.

Orders (800) 252-2828 • Technical Hotline (877) 242-2423 • www.daw-usa.com




PROPER FITTING

Optimum results can only be obtained when DAWFlo Sock and Sheath are fitted properly to provide gradient compression. The WIDTH of the unit MUST BE SMALL ENOUGH TO EFFECTIVELY USE THE ELASTIC QUALITIES IN THE WIDTH (not length) OF THE SOCK AND SHEATH. If fitted incorrectly, the patient may experience some good, but not excellent results.



CORRECT FIT

DAWFlo sock and Sheath are snug, smooth and wrinkle free. WIDTH SIZE IS CORRECT AND UNITS WILL NOT SLIP DOWN. Sufficient compression is applied.



INCORRECT FIT

Unit may be wrinkle free and appear snug, but width is not small enough. DAWFlo Sheath and Sock will not provide sufficient compression.

If you question whether your patient is getting the best fit, please:


- Follow instructions on next page to obtain necessary measurements.
- Call us with the measurements, toll free at (800) 252-2828.



**DAW[®]
INDUSTRIES**


HOW TO MEASURE FOR PROPER DAWFLO® FITTING

ARM




- A. Measure the circumference on the limb at the level of the proximal socket brim. =
- B. Measure the circumference 2 inches above the distal end of the residual limb. =
- C. Measure the length from the level where the proximal circumference measurement was taken down to the distal end of the residual limb. =

A.K.




- A. Measure the circumference on the limb from the top of the lateral wall to the top of the medial wall (this is different from the ischial circumference). =
- B. Measure the circumference at the ischial (crotch) level. =
- C. Measure the circumference 2 inches above the distal end of the residual limb. =
- D. Measure the length from the top of the lateral wall. =
- E. Measure distance between top of Lateral wall to ischial level (crotch). =

B.K.



- A. Measure the circumference on the limb at the level of the top of the suspension strap if you are using one. =
- B. Measure at the level of the top of socket brim. =
- C. Measure the circumference 2 inches above the distal end of the residual limb. =
- D. Measure the length from top of socket down to the distal end of the residual limb. =
- E. In case you wear a Suspension Strap: measure distance between top of Suspension Strap and top of socket. =

SYMES



- A. Measure the circumference on the limb at the level of the top of the suspension strap if you are using one. =
- B. Measure at the level of the top of socket brim. =
- C. Measure the circumference where ankle is narrow. =
- D. Measure circumference at largest distal end. =
- E. Measure the length from top of socket down to the distal end of the residual limb. =
- F. In case you wear a Suspension Strap: measure distance between top of Suspension Strap and top of socket. =

Please call us at 800-252-2828 with measurements.



**DAW®
INDUSTRIES**

© Copyright 2015, DAW Industries, San Diego, CA. All Rights Reserved.

Orders (800) 252-2828 • Technical Hotline (877) 242-2423 • www.daw-usa.com





DAWFLO™

Here is an example of an order with all the information needed:

QUANTITY	DESCRIPTION	SIZE	DESIGN
12	3-Ply DAWFLO Ensemble	1-12	Yellow
6	5-Ply DAWFLO Ensemble	1-14	Green
3	DAWFLO Sheaths	2-16	Yellow
2	5-Ply DAWFLO Sock	3-18	

For AK no need to specify design.

1

QUANTITY

Give the quantity for each item and size you wish to order. If your total order includes 48 packs or more you are entitled to the discount price in the "48 or more items" column.

2

DESCRIPTION *(If not specified, 5-Ply combination packs will be shipped).*

Tell us if you wish the pack to include:

DAW SHEATH ONLY.....	includes one DAW sheath with ventilated mesh
OR A 3-PLY SOCK ONLY.....	includes one fleeced 3-ply DAW sock. *
OR A 5-PLY SOCK ONLY.....	includes one fleeced 5-ply DAW sock. *
OR A 6-PLY SOCK ONLY.....	includes one fleeced 6-ply DAW sock. *
OR A 3-PLY DAW ENSEMBLE.....	includes one DAW sheath and one fleeced 3-ply* DAW sock in the same size as well as a washing bag.
OR A 5-PLY DAW ENSEMBLE.....	includes one DAW sheath and one fleeced 5-ply* DAW sock in the same size as well as a washing bag.
OR A 6-PLY DAW ENSEMBLE.....	includes one DAW sheath and one fleeced 6-ply* DAW sock in the same size as well as a washing bag.

*also available: 3-ply unfleeced socks (thinner than a regular 3-ply sock they correspond in thickness to a 2-ply sock)

3

SIZE

Give: Either the American Size

#1-12 or Top 7in - Toe 4in - 12in long

#13-18 or Top 9in - Toe 6in - 18in long

The Metric Size as follows:

= F2-3Ø

EXAMPLE: F4-F5

Caution:

DAW socks can be machine washed and dried, they do not shrink, they keep elasticity and softness—Some patients who are used to using other manufacturers' socks might be wearing a size larger than needed considering the shrinkage. With DAWFLO, to enhance venous return, you want a snug fit from the first fitting

—For example you might want to use a size 1 Instead of a size 2 for such a patient.



**DAW®
INDUSTRIES**


HOW TO ORDER

4


DESIGN SELECTION FOR DAWFLOW SHEATH & SOCK

- For AK as well as BK with thigh corset, no need to specify the design (it will always be supplied in the yellow design).
- For BK amputees in American size 1 and size 2 you have the choice between the yellow design and the green design:

Up to 18 inches (45 cm) long size 1 and size 2, the socks and sheaths in the yellow design have the same shape as other manufacturers' prosthetic socks, in the green design they have a more tapered shape and are for BK with a pointed end of residual limb



20 inches (50 cm) or longer, the green design offers a bulbous distal end:



DAWFLO SHEATH AND SOCK WITH GARTER

YELLOW DESIGN with V-Band Garter (in yellow pack with a blue stripe)See Yellow Design with V-Band Garter
YELLOW DESIGN with Garter (in yellow pack with a blue stripe)See Yellow Design schedule (non-inventory item*)
GREEN DESIGN with Garter (in green packs with a blue stripe)Please Call Customer Service (non-inventory item*)

- CUSTOM ORDERS -
SEND TRACING

35% ADDITIONAL
MINIMUM QUANTITY: 6

*Custom Order and Non-inventory Items
Are Not Returnable



**DAW®
INDUSTRIES**

© Copyright 2015, DAW Industries, San Diego, CA. All Rights Reserved.

Orders (800) 252-2828 • Technical Hotline (877) 242-2423 • www.daw-usa.com





SIZING CHART

Worldly acclaimed as the finest sock / sheath system, it provides the ultimate care of the residual limb by helping promote venous return. A technical breakthrough together with old-fashioned workmanship in handsewing brings a matchless level of quality. The DAWFLO sock / sheath ensemble is one of the greatest advancements in prosthetic care - it is not "just another sock".



DAW[®]
INDUSTRIES



DAWFLO Sock Only

FOO	DAWFLO System™			DAWFLO Sock Only		
	Qty	3-PLY Ensemble	Qty	5-PLY Ensemble	Stock#	Stock#
Yellow Design						
LENGTH	TOP WIDTH LYING FLAT	TOE WIDTH LYING FLAT	STOCK#	STOCK#	STOCK#	STOCK#
8" (20cm)	3.50"	2.50"	300C51 F00-20	500C51 F00-20	10-01 F00-20	50-51 F00-20
10" (25cm)	3.50"	2.50"	300C51 F00-25	500C51 F00-25	10-01 F00-25	50-51 F00-25
12" (30cm)	3.50"	2.25"	300C51 F00-30	500C51 F00-30	10-01 F00-30	50-51 F00-30
14" (35cm)	3.50"	2.25"	300C51 F00-35	500C51 F00-35	10-01 F00-35	50-51 F00-35
16" (40cm)	3.50"	2.25"	300C51 F00-40	500C51 F00-40	10-01 F00-40	50-51 F00-40
18" (45cm)	3.50"	2.25"	300C51 F00-45	500C51 F00-45	10-01 F00-45	50-51 F00-45
20" (50cm)	3.50"	2.25"	300C51 F00-50	500C51 F00-50	10-01 F00-50	50-51 F00-50

Qty = Also available in 100% Cotton. Stock numbers change to: 30C90, 53-90, etc. Price is 10% less.



DAWFLO SIZING CHART

SIZE A		DAWFLO System™		DAWFLO Sock Only	
Yellow Design		Gtt 3-PLY Ensemble	Gtt 5-PLY Ensemble	6-PLY Ensemble	Gtt 3-PLY Sock Only
LENGTH	TOP WIDTH LYING FLAT	STOCK#	STOCK#	STOCK#	STOCK#
8" (20cm)	4.00"	2.75"	30C51 F0-20	50C51 F0-20	30-01 F0-20
10" (25cm)	4.00"	2.75"	30C51 F0-25	50C51 F0-25	30-01 F0-25
12" (30cm)	4.00"	2.75"	30C51 F0-30	50C51 F0-30	30-01 F0-30
14" (35cm)	4.00"	2.75"	30C51 F0-35	50C51 F0-35	30-01 F0-35
16" (40cm)	4.00"	2.50"	30C51 F0-40	50C51 F0-40	30-01 F0-40
18" (45cm)	4.00"	2.50"	30C51 F0-45	50C51 F0-45	30-01 F0-45
20" (50cm)	4.00"	2.50"	30C51 F0-50	50C51 F0-50	30-01 F0-50

Gtt = Also available in 100% Cotton. Stock numbers change to: 30C90, 53-90, etc. Price is 10% less.

SIZE B		DAWFLO System™		DAWFLO Sock Only	
Yellow Design		Gtt 3-PLY Ensemble	Gtt 5-PLY Ensemble	6-PLY Ensemble	Gtt 3-PLY Sock Only
LENGTH	TOP WIDTH LYING FLAT	STOCK#	STOCK#	STOCK#	STOCK#
8" (20cm)	5.00"	3.00"	30C51 F1-20	N/A	10-01 F1-20
10" (25cm)	5.00"	3.00"	30C51 F1-25	50C51 F1-25	10-01 F1-25
12" (30cm)	5.00"	3.00"	30C51 F1-30	50C51 F1-30	10-01 F1-30
14" (35cm)	5.00"	3.00"	30C51 F1-35	50C51 F1-35	10-01 F1-35
16" (40cm)	5.00"	3.00"	30C51 F1-40	50C51 F1-40	10-01 F1-40
18" (45cm)	5.00"	3.00"	30C51 F1-45	50C51 F1-45	10-01 F1-45
20" (50cm)	5.00"	3.00"	30C51 F1-50	50C51 F1-50	10-01 F1-50

Gtt = Also available in 100% Cotton. Stock numbers change to: 30C90, 53-90, etc. Price is 10% less.





SIZING CHART

SIZE Ø Yellow Design (Standard Toe)

DAWFLO System™
Includes: DAWFLO Sheath, DAWFLO Sock, and Washing Bag

DAWFLO Only		DAWFLO Sock Only	
Grt	3-PLY Ensemble	Grt	5-PLY Ensemble
STOCK#	STOCK#	STOCK#	STOCK#
LENGTH	TOP WIDTH LYING FLAT	TOE WIDTH LYING FLAT	
8" (20cm)	5.00"	3.00"	30C51 F1-20
10" (25cm)	5.00"	3.00"	30C51 F1-25
12" (30cm)	5.00"	3.00"	30C51 F1-30
14" (35cm)	5.00"	3.00"	30C51 F1-35
16" (40cm)	5.00"	3.00"	30C51 F1-40
18" (45cm)	5.00"	3.00"	30C51 F1-45
20" (50cm)	5.00"	3.00"	30C51 F1-50
22" (55cm)	*	*	*
24" (60cm)	*	*	*
26" (65cm)	*	*	*
28" (70cm)	*	*	*

* = Also available in 100% Cotton. Stock numbers change to: 30C90, 53-90, etc. Price is 10% less.

SIZE Ø Green Design (Narrow Toe)

DAWFLO System™
Includes: DAWFLO Sheath, DAWFLO Sock, and Washing Bag

DAWFLO Sheath		DAWFLO Sock Only	
Grt	3-PLY Ensemble	Grt	5-PLY Sock Only
STOCK#	STOCK#	STOCK#	STOCK#
LENGTH	TOP WIDTH LYING FLAT	TOE WIDTH LYING FLAT	
8" (20cm)	5.00"	3.25"	35C51 F1-20
10" (25cm)	5.00"	3.25"	35C51 F1-25
12" (30cm)	5.00"	3.25"	35C51 F1-30
14" (35cm)	5.25"	3.00"	35C51 F1-35
16" (40cm)	5.00"	3.00"	35C51 F1-40
18" (45cm)	5.25"	3.00"	35C51 F1-45
20" (50cm)	5.00"	3.00"	35C52 F1-50
22" (55cm)	5.25"	3.00"	35C52 F1-55
24" (60cm)	5.00"	3.00"	35C52 F1-60
26" (65cm)	5.00"	3.00"	35C52 F1-65
28" (70cm)	*	*	*

* = Special Order - Non-Returnable Item.

Grt = Also available in 100% Cotton. Stock numbers change to: 30C90, 53-90, etc. Price is 10% less.

* = Special Order - Non-Returnable Item.





SIZING CHART

SIZE 1		DAWFLO System™		DAWFLO Sheath Only		DAWFLO Sock Only	
	Yellow Design	Grt	3-PLY Ensemble	Grt	5-PLY Ensemble	Grt	6-PLY Ensemble
LENGTH	TOP WIDTH LYING FLAT	STOCK#	STOCK#	STOCK#	STOCK#	STOCK#	STOCK#
8"	(20cm) 6.50"	33C50 F2-20	53C50 F2-20	63C50 F2-20	13-00 F2-20	33-50 F2-20	63-50 F2-20
10"	(25cm) 6.25"	33C50 F2-25	53C50 F2-25	63C50 F2-25	13-00 F2-25	33-50 F2-25	63-50 F2-25
12"	(30cm) 6.50"	33C50 F2-30	53C50 F2-30	63C50 F2-30	13-00 F2-30	33-50 F2-30	63-50 F2-30
14"	(35cm) 6.50"	33C50 F2-35	53C50 F2-35	63C50 F2-35	13-00 F2-35	33-50 F2-35	63-50 F2-35
16"	(40cm) 6.50"	33C50 F2-40	53C50 F2-40	63C50 F2-40	13-00 F2-40	33-50 F2-40	63-50 F2-40
18"	(45cm) 6.50"	33C50 F2-45	53C50 F2-45	63C50 F2-45	13-00 F2-45	33-50 F2-45	63-50 F2-45
20"	(50cm) 6.50"	33C50 F2-50	53C50 F2-50	63C50 F2-50	13-00 F2-50	33-50 F2-50	63-50 F2-50
22"	(55cm) 6.00"	33C51 F2-55	53C51 F2-55	63C51 F2-55	13-01 F2-55	33-51 F2-55	63-51 F2-55
24"	(60cm) 6.00"	33C51 F2-60	53C51 F2-60	63C51 F2-60	13-01 F2-60	33-51 F2-60	63-51 F2-60
26"	(65cm) 6.00"	33C51 F2-65	*	63C51 F2-65	13-01 F2-65	33-51 F2-65	63-51 F2-65
28"	(70cm) *	*	*	*	*	*	*

NOTE: Also available see Section 1.1E p.03

Grt = Special Order – Non-Returnable Item.

* = V-Band Garter available in 100% Cotton. Stock numbers change to: 30C90, 53-90, etc. Price is 10% less.

SIZE 1		DAWFLO System™		DAWFLO Sheath Only		DAWFLO Sock Only	
	Green Design (Narrow Toe)	Grt	3-PLY Ensemble	Grt	5-PLY Ensemble	Grt	6-PLY Ensemble
LENGTH	TOP WIDTH LYING FLAT	STOCK#	STOCK#	STOCK#	STOCK#	STOCK#	STOCK#
8"	(20cm) 6.00"	35C51 F2-20	55C51 F2-20	65C51 F2-20	15-01 F2-20	35-51 F2-20	65-51 F2-20
10"	(25cm) 6.00"	35C51 F2-25	55C51 F2-25	65C51 F2-25	15-01 F2-25	35-51 F2-25	65-51 F2-25
12"	(30cm) 6.00"	35C51 F2-30	55C51 F2-30	65C51 F2-30	15-01 F2-30	35-51 F2-30	65-51 F2-30
14"	(35cm) 6.00"	35C51 F2-35	55C51 F2-35	65C51 F2-35	15-01 F2-35	35-51 F2-35	65-51 F2-35
16"	(40cm) 6.00"	35C51 F2-40	55C51 F2-40	65C51 F2-40	15-01 F2-40	35-51 F2-40	65-51 F2-40
18"	(45cm) 6.00"	35C51 F2-45	55C51 F2-45	65C51 F2-45	15-01 F2-45	35-51 F2-45	65-51 F2-45
20"	(50cm) 6.50"	35C52 F2-50	55C52 F2-50	65C52 F2-50	15-02 F2-50	35-52 F2-50	65-52 F2-50
22"	(55cm) 6.50"	35C52 F2-55	55C52 F2-55	65C52 F2-55	15-02 F2-55	35-52 F2-55	65-52 F2-55
24"	(60cm) 6.00"	35C52 F2-60	55C52 F2-60	65C52 F2-60	15-02 F2-60	35-52 F2-60	65-52 F2-60
26"	(65cm) 6.00"	35C52 F2-65	*	65C52 F2-65	15-02 F2-65	35-52 F2-65	65-52 F2-65
28"	(70cm) *	*	*	*	*	*	*

Grt = Also available in 100% Cotton. Stock numbers change to: 30C90, 53-90, etc. Price is 10% less.

* = Special Order – Non-Returnable Item.



DAWFLO SIZING CHART

SIZE 2		DAWFLO System™		DAWFLO Sheath Only		DAWFLO Sock Only	
Gtt	3-PLY Ensemble	Gtt	5-PLY Ensemble	Gtt	3-PLY Sock Only	Gtt	5-PLY Sock Only
LENGTH	TOP WIDTH LYING FLAT	STOCK#	STOCK#	STOCK#	STOCK#	STOCK#	STOCK#
8" (20cm)	7.50"	4.50"	33C50 F3-20	63C50 F3-20	13-00 F3-20	33-50 F3-20	63-50 F3-20
10" (25cm)	7.50"	4.50"	33C50 F3-25	63C50 F3-25	13-00 F3-25	33-50 F3-25	63-50 F3-25
12" (30cm)	7.50"	4.50"	33C50 F3-30	63C50 F3-30	13-00 F3-30	33-50 F3-30	63-50 F3-30
14" (35cm)	7.50"	4.25"	33C50 F3-35	63C50 F3-35	13-00 F3-35	33-50 F3-35	63-50 F3-35
16" (40cm)	7.50"	4.25"	33C50 F3-40	63C50 F3-40	13-00 F3-40	33-50 F3-40	63-50 F3-40
18" (45cm)	7.50"	4.00"	33C50 F3-45	63C50 F3-45	13-00 F3-45	33-50 F3-45	63-50 F3-45
20" (50cm)	7.50"	3.75"	33C50 F3-50	63C50 F3-50	13-00 F3-50	33-50 F3-50	63-50 F3-50
22" (55cm)	7.50"	3.75"	33C51 F3-55	63C51 F3-55	13-01 F3-55	33-51 F3-55	63-51 F3-55
24" (60cm)	7.25"	3.50"	33C51 F3-60	63C51 F3-60	13-01 F3-60	33-51 F3-60	63-51 F3-60
26" (65cm)	7.25"	3.50"	33C51 F3-65	63C51 F3-65	13-01 F3-65	33-51 F3-65	63-51 F3-65
28" (70cm)	*	*	*	*	*	*	*

SIZE 2		DAWFLO System™		DAWFLO Sheath Only		DAWFLO Sock Only	
Gtt	3-PLY Ensemble	Gtt	5-PLY Ensemble	Gtt	3-PLY Sock Only	Gtt	5-PLY Sock Only
LENGTH	TOP WIDTH LYING FLAT	STOCK#	STOCK#	STOCK#	STOCK#	STOCK#	STOCK#
8" (20cm)	7.00"	4.00"	35C51 F3-20	65C51 F3-20	15-01 F3-20	35-51 F3-20	65-51 F3-20
10" (25cm)	7.00"	4.00"	35C51 F3-25	65C51 F3-25	15-01 F3-25	35-51 F3-25	65-51 F3-25
12" (30cm)	7.00"	4.00"	35C51 F3-30	65C51 F3-30	15-01 F3-30	35-51 F3-30	65-51 F3-30
14" (35cm)	6.75"	4.00"	35C51 F3-35	65C51 F3-35	15-01 F3-35	35-51 F3-35	65-51 F3-35
16" (40cm)	6.75"	3.50"	35C51 F3-40	65C51 F3-40	15-01 F3-40	35-51 F3-40	65-51 F3-40
18" (45cm)	7.00"	3.50"	35C51 F3-45	65C51 F3-45	15-01 F3-45	35-51 F3-45	65-51 F3-45
20" (50cm)	6.75"	3.50"	35C52 F3-50	65C52 F3-50	15-02 F3-50	35-52 F3-50	65-52 F3-50
22" (55cm)	6.75"	3.50"	35C52 F3-55	65C52 F3-55	15-02 F3-55	35-52 F3-55	65-52 F3-55
24" (60cm)	6.75"	3.50"	35C52 F3-60	65C52 F3-60	15-02 F3-60	35-52 F3-60	65-52 F3-60
26" (65cm)	6.50"	*	*	*	*	*	*
28" (70cm)	*	*	*	*	*	*	*

V = V-Band Garter available see Section 1.1E p.03

U = Also available in Knee Disarticulation. Stock # changes to 33C53 or 53C53 or 13-03

Gtt = Also available in 100% Cotton. Stock numbers change to: 30C90, 53-90, etc. Price is 10% less.

* = Special Order - Non-Returnable Item.

V = V-Band Garter available see Section 1.1E p.03
 Gtt = Also available in 100% Cotton. Stock numbers change to: 30C90, 53-90, etc. Price is 10% less.
 * = Special Order - Non-Returnable Item.





SIZING CHART

DAWFLO System™

Includes: DAWFLO Sheath, DAWFLO Sock, and Washing Bag

SIZE 3

Yellow Design

DAWFLO System™		DAWFLO Sock Only						
DAWFLO Sheath		Only	Ctt	3-PLY Sock Only	Ctt	5-PLY Sock Only	Ctt	6-PLY Sock Only
		STOCK#	STOCK#	STOCK#	STOCK#	STOCK#	STOCK#	STOCK#
Ø Ø	8" (20cm)	33C50 F4-20	53C50 F4-20	13-00 F4-20	33-50 F4-20	53-50 F4-20	63-50 F4-20	
Ø Ø	10" (25cm)	33C50 F4-25	53C50 F4-25	13-00 F4-25	33-50 F4-25	53-50 F4-25	63-50 F4-25	
Ø Ø	12" (30cm)	33C50 F4-30	53C50 F4-30	13-00 F4-30	33-50 F4-30	53-50 F4-30	63-50 F4-30	
Ø Ø	14" (35cm)	33C50 F4-35	53C50 F4-35	13-00 F4-35	33-50 F4-35	53-50 F4-35	63-50 F4-35	
U U	16" (40cm)	33C50 F4-40	53C50 F4-40	13-00 F4-40	33-50 F4-40	53-50 F4-40	63-50 F4-40	
U U	18" (45cm)	33C50 F4-45	53C50 F4-45	13-00 F4-45	33-50 F4-45	53-50 F4-45	63-50 F4-45	
U U	20" (50cm)	33C50 F4-50	53C50 F4-50	13-00 F4-50	33-50 F4-50	53-50 F4-50	63-50 F4-50	
U U	22" (55cm)	33C51 F4-55	53C51 F4-55	13-01 F4-55	33-51 F4-55	53-51 F4-55	63-51 F4-55	
U U	24" (60cm)	33C51 F4-60	53C51 F4-60	13-01 F4-60	33-51 F4-60	53-51 F4-60	63-51 F4-60	
U U	26" (65cm)	33C51 F4-65	53C51 F4-65	13-01 F4-65	33-51 F4-65	53-51 F4-65	63-51 F4-65	
U U	28" (70cm)	*	*	*	*	*	*	*

Ø = V-Band Garter available see Section 1.1E p.03

U = Also available in Knee Disarticulation. Stock # changes to 33C53 or 53C53 or 13-03

Ctt = Also available in 100% Cotton. Stock numbers change to: 30C90, 53-90, etc. Price is 10% less.

* = Special Order - Non-Returnable Item.

DAWFLO System™

Includes: DAWFLO Sheath, DAWFLO Sock, and Washing Bag

SIZE 4

Yellow Design

DAWFLO System™		DAWFLO Sock Only						
DAWFLO Sheath		Only	Ctt	3-PLY Sock Only	Ctt	5-PLY Sock Only	Ctt	6-PLY Sock Only
		STOCK#	STOCK#	STOCK#	STOCK#	STOCK#	STOCK#	STOCK#
Ø Ø	8" (20cm)	33C50 F5-20	53C50 F5-20	13-00 F5-20	33-50 F5-20	53-50 F5-20	63-50 F5-20	
Ø Ø	10" (25cm)	33C50 F5-25	53C50 F5-25	13-00 F5-25	33-50 F5-25	53-50 F5-25	63-50 F5-25	
Ø Ø	12" (30cm)	33C50 F5-30	53C50 F5-30	13-00 F5-30	33-50 F5-30	53-50 F5-30	63-50 F5-30	
Ø Ø	14" (35cm)	33C50 F5-35	53C50 F5-35	13-00 F5-35	33-50 F5-35	53-50 F5-35	63-50 F5-35	
U U	16" (40cm)	33C50 F5-40	53C50 F5-40	13-00 F5-40	33-50 F5-40	53-50 F5-40	63-50 F5-40	
U U	18" (45cm)	33C50 F5-45	53C50 F5-45	13-00 F5-45	33-50 F5-45	53-50 F5-45	63-50 F5-45	
U U	20" (50cm)	33C50 F5-50	53C50 F5-50	13-00 F5-50	33-50 F5-50	53-50 F5-50	63-50 F5-50	
U U	22" (55cm)	33C51 F5-55	53C51 F5-55	13-01 F5-55	33-51 F5-55	53-51 F5-55	63-51 F5-55	
U U	24" (60cm)	33C51 F5-60	53C51 F5-60	13-01 F5-60	33-51 F5-60	53-51 F5-60	63-51 F5-60	
U U	26" (65cm)	33C51 F5-65	53C51 F5-65	13-01 F5-65	33-51 F5-65	53-51 F5-65	63-51 F5-65	
U U	28" (70cm)	*	*	*	*	*	*	*

Ø = V-Band Garter available see Section 1.1E p.03

U = Also available in Knee Disarticulation. Stock # changes to 33C53 or 53C53 or 13-03

Ctt = Also available in 100% Cotton. Stock numbers change to: 30C90, 53-90, etc. Price is 10% less.

* = Special Order – Non-Returnable Item.





D A W F L O S I Z I N G C H A R T

Size 1 V-BAND GARTER Yellow Design 7" TOP-4" TOE-LYING FLAT

DAWFLO System™				DAWFLO Sock Only			
Qty:	3-PLY Ensemble	Qty:	5-PLY Ensemble	Qty:	DAWFLO Sheath Only	Qty:	5-PLY Sock Only
LENGTH	TOP WIDTH LYING FLAT	TOE WIDTH LYING FLAT	STOCK#	STOCK#	STOCK#	STOCK#	STOCK#
8" (20cm)	6.50"	4.00"	34C54 F2-20	54C54 F2-20	* 14-04 F2-20	34-54 F2-20	54-54 F2-20
10" (25cm)	6.25"	4.00"	34C54 F2-25	54C54 F2-25	* 14-04 F2-25	34-54 F2-25	54-54 F2-25
12" (30cm)	6.50"	4.50"	34C54 F2-30	54C54 F2-30	* 14-04 F2-30	34-54 F2-30	54-54 F2-30
14" (35cm)	6.50"	4.00"	34C54 F2-35	54C54 F2-35	* 14-04 F2-35	34-54 F2-35	54-54 F2-35
16" (40cm)	6.50"	3.25"	34C54 F2-40	54C54 F2-40	* 14-04 F2-40	34-54 F2-40	54-54 F2-40
18" (45cm)	6.50"	3.25"	34C54 F2-45	54C54 F2-45	* 14-04 F2-45	34-54 F2-45	54-54 F2-45
20" (50cm)	6.50"	3.25"	34C54 F2-50	54C54 F2-50	* 14-04 F2-50	34-54 F2-50	54-54 F2-50

Qty = Also available in 100% Cotton. Stock numbers change to: 30C90, 53-90, etc. Price is 10% less.
 * = Special Order – Non-Returnable Item.

Size 2 V-BAND GARTER Yellow Design 8" TOP-5" TOE-LYING FLAT

DAWFLO System™				DAWFLO Sock Only			
Qty:	3-PLY Ensemble	Qty:	5-PLY Ensemble	Qty:	DAWFLO Sheath Only	Qty:	5-PLY Sock Only
LENGTH	TOP WIDTH LYING FLAT	TOE WIDTH LYING FLAT	STOCK#	STOCK#	STOCK#	STOCK#	STOCK#
8" (20cm)	7.50"	4.50"	34C54 F3-20	54C54 F3-20	* 14-04 F3-20	34-54 F3-20	54-54 F3-20
10" (25cm)	7.50"	4.50"	34C54 F3-25	54C54 F3-25	* 14-04 F3-25	34-54 F3-25	54-54 F3-25
12" (30cm)	7.50"	4.50"	34C54 F3-30	54C54 F3-30	* 14-04 F3-30	34-54 F3-30	54-54 F3-30
14" (35cm)	7.50"	4.25"	34C54 F3-35	54C54 F3-35	* 14-04 F3-35	34-54 F3-35	54-54 F3-35
16" (40cm)	7.50"	4.25"	34C54 F3-40	54C54 F3-40	* 14-04 F3-40	34-54 F3-40	54-54 F3-40
18" (45cm)	7.50"	4.00"	34C54 F3-45	54C54 F3-45	* 14-04 F3-45	34-54 F3-45	54-54 F3-45
20" (50cm)	7.50"	3.75"	34C54 F3-50	54C54 F3-50	* 14-04 F3-50	34-54 F3-50	54-54 F3-50

Qty = Also available in 100% Cotton. Stock numbers change to: 30C90, 53-90, etc. Price is 10% less.
 * = Special Order – Non-Returnable Item.





SIZING CHART

Size 3 V-BAND GARTER Yellow Design 9" TOP-6" TOE-LYING FLAT

DAWFLO System™

Includes: DAWFLO Sheath, DAWFLO Sock, and Washing Bag

Length	Top Width	Toe Width	Lying Flat	DAWFLO Sheath Only			DAWFLO Sock Only		
				Cotton	3-Ply Ensemble	Cotton	5-Ply Ensemble	Cotton	6-Ply Ensemble
			Stock#		Stock#		Stock#		Stock#
8" (20cm)	9.00"	5.00"	34C54 F4-20	*	54C54 F4-20	*	14-04 F4-20	34-54 F4-20	*
10" (25cm)	9.00"	5.00"	34C54 F4-25	*	54C54 F4-25	*	14-04 F4-25	34-54 F4-25	*
12" (30cm)	9.00"	5.00"	34C54 F4-30	*	54C54 F4-30	*	14-04 F4-30	34-54 F4-30	*
14" (35cm)	9.25"	5.00"	34C54 F4-35	*	54C54 F4-35	*	14-04 F4-35	34-54 F4-35	*
16" (40cm)	9.25"	5.00"	34C54 F4-40	*	54C54 F4-40	*	14-04 F4-40	34-54 F4-40	*
18" (45cm)	9.25"	5.00"	34C54 F4-45	*	54C54 F4-45	*	14-04 F4-45	34-54 F4-45	*
20" (50cm)	9.25"	4.00"	34C54 F4-50	*	54C54 F4-50	*	14-04 F4-50	34-54 F4-50	*

Cotton = Also available in 100% Cotton. Stock numbers change to: 30C90, 53-90, etc. Price is 10% less.
 * = Special Order – Non-Returnable Item.

Size 3 V-BAND GARTER Yellow Design 10" TOP-5" TOE-LYING FLAT

DAWFLO System™

Includes: DAWFLO Sheath, DAWFLO Sock, and Washing Bag

Length	Top Width	Toe Width	Lying Flat	DAWFLO Sheath Only			DAWFLO Sock Only		
				Cotton	3-Ply Ensemble	Cotton	5-Ply Ensemble	Cotton	6-Ply Ensemble
			Stock#		Stock#		Stock#		Stock#
8" (20cm)	10.25"	5.75"	34C54 F5-20	*	54C54 F5-20	*	14-04 F5-20	34-54 F5-20	*
10" (25cm)	10.25"	5.75"	34C54 F5-25	*	54C54 F5-25	*	14-04 F5-25	34-54 F5-25	*
12" (30cm)	10.25"	5.75"	34C54 F5-30	*	54C54 F5-30	*	14-04 F5-30	34-54 F5-30	*
14" (35cm)	10.25"	5.75"	34C54 F5-35	*	54C54 F5-35	*	14-04 F5-35	34-54 F5-35	*
16" (40cm)	10.25"	5.75"	34C54 F5-40	*	54C54 F5-40	*	14-04 F5-40	34-54 F5-40	*
18" (45cm)	10.25"	5.75"	34C54 F5-45	*	54C54 F5-45	*	14-04 F5-45	34-54 F5-45	*
20" (50cm)	10.25"	5.25"	34C54 F5-50	*	54C54 F5-50	*	14-04 F5-50	34-54 F5-50	*

Cotton = Also available in 100% Cotton. Stock numbers change to: 30C90, 53-90, etc. Price is 10% less.
 * = Special Order – Non-Returnable Item.



DAW®
INDUSTRIES

