

# DAWFLO<sup>TM</sup> SHEATH



## CLAIM TO FAME

- Helps venous return, capillary stimulation when properly fitted.
- Reduces or eliminates irritation, friction and bruising.
- Moisture barrier, substantially decreases perspiration.
- Helps keep skin dry.
- Circulation enhancement through gradient pressure.
- Outlives regular Prosthetic Sheaths 10 to 1.

Worldly acclaimed as the finest Prosthetic Sheath system, it provides the ultimate care of the residual limb by helping promote venous return. A technical breakthrough together with old-fashioned workmanship in handsewing brings a matchless level of quality. The **DAWFLO** Sheath is one of the greatest advancements in prosthetic care - it is not "just another prosthetic sheath".

The **DAWFLO** Sheath has a soft surface next to the skin made of the finest sheer material. Its outside surface grips the **DAWFLO** Sock, preventing movement between Sheath and Sock. The perspiration passes through to the **DAWFLO** Sock thanks to the **DAWFLO** Sheath ventilated mesh.

When properly fitted, it supports the skin and muscle tissue of the residual limb. The **DAWFLO** Sheath is highly effective in helping to create a perfect, total fit of the prosthesis to the limb. The special design allows for easy circulation of the air, thus evaporating the perspiration and creating a moisture barrier screen which helps keep the skin dry.

## HOW TO ORDER

Please refer to the chart in **Section 1.1F, p.02**



**DAW<sup>®</sup>  
INDUSTRIES**

