AK ADULT KNEE



Stealth MPK 3.0

#: TGK-5PS10SLKMTX

The technology of the Stealth MPK works like an ultra fast 3D gyroscope which instantly detects the knee positioning in space and reacts by providing through its microprocessor the optimum extension flexion resistance while walking down an incline or stairs. Select on the DAW website a predetermined set of parameters from a menu according to patient weight, level of activity, muscle strength, and preferred activities. Then tweak the parameters to reach that perfect natural gait that everyone hopes for but cannot quite reach.

Benefits:

- ✓ MPK with pinpoint programmable gait accuracy
- ✔ Advanced Pneumatics with extremely high flexion resistance
- ✓ 5-bar stability, fully adjustable stance flexion
- ✓ Full range adjustable geometric stance control
- ✓ Trouble free low maintenance
- ✓ Proven durability and dependability
- ✓ Fully programmable centrode (PC)
- ✓ Recordable gait data for seamless Bluetooth transfer to another Stealth MPK



Specifications:

Function Level	K3/K4		
Indicates	Transfemoral or Hip Disarticulation		
Item Weight	38.6 oz 1095 g		
Weight Limit	275 lb 145 kg		
Pylon diameter	30 mm		
Maximum Flexion Angle	140°		
Water Exposure	Splash Resistant		
Warranty	1-year Optional 3-year Available		


Connectors and Accessories:

Proximal Attachments	Next Page	
ePAD™ Alignment System		p.01
Pylon: (ø 30mm)		p.01
Ankles	Section 6A	p.01
Feet	Section 7A	p.01
Foam Covers: – 1-part DAWFlex™ Cover	Section 9.1A	p.01
– Discontinuous 2-part		p.04
Skin and other finishings		p.01



Orders (800) 252-2828 • Technical Support (800) 242-8669 • www.daw-usa.com

AK ADULT KNEES PROXIMAL ATTACHMENT SELECTOR CHART





Orders (800) 252-2828 • Technical Support (800) 242-8669 • www.daw-usa.com