

HEIGHT ADJUSTABLE PYLONS



Stock Number							
TTT - A6/9	AA	Ø 30mm	Ø 30mm	8 7/8 oz (251 g)	6 1/2 to 8 1/2 in (16.5 to 21.5 cm)	DO NOT DISPENSE Clinical Evaluation Tool Only	



Stock Number							
TTT - A6/9	AA	Ø 30mm	Ø 30mm	11 5/8 oz (330 g)	8 1/2 to 12 1/2 in (21.5 to 31.8 cm)	DO NOT DISPENSE Clinical Evaluation Tool Only	

SYMBOL LEGEND™

Material		ATTACHMENT METHOD		Dimensions		Maximum Patient Weight	
		= Top	= Bottom				
A = Aluminum Alloy	P = Plastic	= 1 Center Hole	= Clamp	= Male Pyl.	= Adhesive	= Bolt Size	= 80 lb (36 kg)
AA = Aerospace Aluminum	S = Stainless	= Euro 4-Hole	= Female Pyl.	= Lamination	= Applied Ht. or Thickness	= 120 lb (54 kg)	= 220 lb (100 kg)
F = Foam	T = Titanium	= 30mm Pylon	= Thermoplastic	= Ht.	= Wt.	= Dia.	= 243 lb (110 kg)
G = Graph-Lite™	W = Wood	= Slide Adjust		= Ht.	= Wt.	= Dia.	= 165 lb (75 kg)
				= Ht.	= Wt.	= Dia.	= 250 lb (113 kg)
				= Ht.	= Wt.	= Dia.	= 200 lb (90 kg)



DAW
INDUSTRIES

