

## KNEE ROTATOR

TKR-Ø1

**Claim to fame**

- Automatic Lock at neutral position
- Can be used with any DAW knee unit
- Perfect match for the Pyramid System
- Fits straight into any DAW Type "T" Knee\* or any standard Euro 4-hole pattern

Crossing one's legs is not just a luxury, it's a necessity. Patients will tell you that this feature is convenient for changing shoes, getting out of the car, sitting on the floor and numerous other activities.

Why not give your patients the opportunity to experience these "luxuries" we take for granted? All it takes is this Knee Rotator and the push of a button.

Stock Number							
TKR-Ø1	AA	4 THREADED (6 mm)	1 7/16 THREADED (6 mm)	6 5/8 oz (188 g)	3/4 in (1.9 cm)		R

\* Type "T" Knees = DAW Adult Knees which stock # starts with a "T" (except TK-40C & TK-1C1)

## SYMBOL LEGEND™

Material		ATTACHMENT METHOD		Dimensions		Maximum Patient Weight	
A = Aluminum Alloy AA = Aerospace Aluminum F = Foam G = Graph-Lite™	P = Plastic S = Stainless T = Titanium W = Wood	⊙ = 1 Center Hole ⊕ = Euro 4-Hole	⊙ = Clamp ⊙ = Male Pyr. ⊙ = Female Pyr. ⊙ = Thermoplastic	L = Lamination A = Adhesive B = Bolt S = Size	T = Ht. Wt. = Wt. D = Dia. A = Applied Ht. or Thickness C = Circ.	⊙ = 80 lb (36 kg) ⊙ = 120 lb (54 kg) ⊙ = 165 lb (75 kg) ⊙ = 200 lb (90 kg)	⊙ = 220 lb (100 kg) ⊙ = 243 lb (110 kg) ⊙ = 250 lb (113 kg) ⊙ = 250 lb (113 kg) high activity



**DAW**  
INDUSTRIES