

MULTI-AXIS FEET

1



Description	Stock Number	Function	Toes	Special Feature(s)					Size R or L	Heel Height			
Foot	TFF-Ø1XXR/L	SACH	Fits DAWtoeSKIN™ Style "A"	Ultra Light	TGA-ØM	TGP-ØM	incl.		22-29	3/4 in (19 mm)	15 ½ oz* (439 oz)		
Set of 2 Bumpers	TFB-N1 or TFB-H1	Normal or High Density						* Foot weight based on size 25 foot XX = SIZE					

Must be used with:

Graph-Lite™ Multi-Axis Ankle
Stock # TGA-ØM


or

Graph-Lite™ Multi-Axis Ankle
Stock # TGP-ØMFC



For additional information, please refer to Section 6A, p.03.

2



Can be used with:

Multi-Axis Ankles
Stock # TAJ-CLNUTxxx
or TAJ-CLØØ2xxx

or

Multi-Axis Ankles
Stock # TAJ-MPYRNUTxxx
or TAJ-MPYRØØ2xxx

* Patients up to 130 lbs (60 kg) xxx = Ø60
" " " 175 lbs (80 kg) xxx = Ø80
" " " 220 lbs (100 kg) xxx = 1Ø0

For additional information, please refer to Section 6A, p.03.

SYMBOL LEGEND™

Material		Alignable	ATTACHMENT METHOD						Dimensions	Maximum Patient Weight			
			= Top			= Bottom							
= Aluminum Alloy	= Plastic	= Angular Adj.	= 1 Center Hole	= Clamp	= Male Pyl.	= Lamination	= Adhesive	= Ht.	= Wt.	= Dib.	= 80 lb (36 kg)	= 220 lb (100 kg)	
= Aerospace Aluminum	= Stainless	= Rotation	= Euro 4-Hole	= 30mm Pylon	= Female Pyr.	= Thermoplastic	= Bolt Size	= Applied Ht. or Thickness	= Circ.	= 120 lb (54 kg)	= 243 lb (110 kg)	= 165 lb (75 kg)	= 250 lb (113 kg)
= Foam	= Titanium	= Slide Adjust									= 200 lb (90 kg)	= high activity	
= Graph-Lite™	= Wood												



DAW INDUSTRIES

