

# COSMETIC STOCKINGS

## COSMETIC STOCKINGS TO PROTECT AND BEAUTIFY



DAW  
Sheer Stocking

DAW  
Ribbed Stocking

### Claim to fame

- Easy to apply
- Very inexpensive
- 3 lengths AK and BK fit all
- Tubular design fits any size foot



The DAW Ribbed Stocking substantially prolongs the life of any Foam Cover. Apply the ribbed stocking(s) first. Then apply sheer stocking, which helps to achieve a more natural looking prosthesis. Since these are tube stockings, please take into account the foot length when measuring your patient.

### HOW TO ORDER

**TCS-SA<sub>xx</sub>**  
**TCS-SB<sub>xx</sub>**

**TCS-RA<sub>xx</sub>**  
**TCS-RB<sub>xx</sub>**

#### Lengths available

**A = AK**..... xx = 35, 38, and 41 in (89, 95, and 104 cm)

**B = BK**..... xx = 18, 21 and 24 in (45, 53, and 60 cm)

xx = Length



**DAW<sup>®</sup>**  
**INDUSTRIES**

© Copyright 2015, DAW Industries, San Diego, CA. All Rights Reserved.

Orders (800) 252-2828 • Technical Hotline (877) 242-2423 • [www.daw-usa.com](http://www.daw-usa.com)